

SETTING UP YOUR TABLE




1. Lay Table on it's Side and Unfold.



2. Open Up the Legs (The Lower Legs in The Original Position When I Start Unfolding).



3. Straighten Each Brace.



4. Lift the Table Up Until All Legs are in the floor.



5-1. Insert the bars of headrest into the heels, put the pillow on the headrest.




5-2. Insert the bars of side support.

LOSE AND TRANSPORT TABLE



- 1. Lay table on its side and fold at an angle.
- 2. Completely fold the legs flat under the table.
- 3. Close the table halves and use the latch to lock the table.


3-SECTION MASSAGE: TABLE



Headrest, Armrest(2PCs), Arm support and Carry Bag are optional part.

PRECAUTIONARY INSTRUCTIONS

1. Fold/unfold the massage table by 2 people.
2. Check & Tighten all the screws before use.
3. Always place the product on flat & stable surface. Make sure all the legs are at the same level to avoid instability, which might cause injury.
4. Maintain the table periodically. Make sure all screws are secured and the cables were not deteriorated or frayed.
5. Do not overload the table beyond the recommended weight limit. This could result in the serious injury to the user. Do not put your fingers into the gap of the table nor the buckle at the bottom of the table.
6. Do not use alcohol-based cleaners when cleaning the PVC cover. This could cause damage to the material. Wipe clean with a damp cloth only.



Part#	Description	Quantity
1.	Leg knob	4
2.	Long cable	2
3.	Leg cushion	4
4.	Hinge	3
5.	Handle	1
6.	PVC cover	3
7.	Bracket	2
8.	Brace	2
9.	Support	4

OVERVIEW & PACKING LIST